GO-TAPER UNIVERSAL

1. Specifications

Model: A0, A1, A2, B1, B2, B3, B4, B5
 Length : 19mm, 21mm, 25mm, 31mm
 Package: Plastic box- 6pcs/box

2. Raw materials

NiTi alloy

3. Function

This is the dental instrument with columnar or tapering working part and it possesses circular form in cross cut section and flat end. GO-TAPER UNIVERSAL files are used aiming at shaping and cleaning the root canal by mechanically powered GO-TAPER UNIVERSAL system.

4. Recommend Speed and Torque

Model	Length (mm)	Taper	Speed (RPM)	Torque
A0	19	4%	300	5
A1	21, 25, 31	2%	300	5
A2&B1	21, 25, 31	A2 : 4% B1 : 7%	300	2
B2	21, 25, 31	B2 : 8%	300	3
B3		B3 : 9%		
B4		B4 : 6%		
B5		B5 : 5%		

5. Motion

Continuous rotation

6. Precautions for use

Safety equipment (gloves, glasses, rubber dam) is required for patients and self-protection.

These instruments can only be used after sterilization, please using these files under medical environment.

Please use GO-TAPER UNIVERSAL system in an appropriate case and in a correct way. Remove debris out frequently from screwing-in problems

Normally shaping the canal with EDTA for entering the canal smoothly and irrigate sufficiently after removing debris created by shaping the root canal.

If you come across a situation that it's difficult to enter the deeper part of the canal, please stop and check both morphology of the canal and whether there are signs of deformation of files or not, then confirm the orifice is expanded wide enough for subsequent use of shaping files.

7. Selecting the correct GO-TAPER UNIVERSAL

A0- Opener

A1, A2 work to smooth canal, shaping canal

B1, B2, B3, B4, B5 are finishing files.

Normally B2 is enough for most of case.

If canal is bigger, then use B3, B4, B5.

8. Sequence

- X-ray film to estimate the working length of canals
- Create straightline access to canal orifice
- In the presence with EDTA, scout the coronal two-thirds with SS hand K-files #10/#15, then gently work K-files until a smooth glide path is confirmed
- Rinsing the canal with NaOCl enough, use A1 and follow the glide path until light resistance, then laterally "brush" and cut dentin on the outstroke
- Shaping with A1 as described until the depth of the #15 hand file is reached
- Use A2, exactly as described for A1, until the depth of the #15 hand file is reached
- In the presence with EDTA, scout the apical one-third with SS hand K-files #10/#15, and gently work them until they are loose at length

- Establish working length, confirm patency and verify a smooth glide path in the apical one-third
- Use A1, as described, until working length is reached
- Use A2, as described, until working length is reached
- Confirm working length
- Use B1 to shape the apical area
- Gauge the shaping canals with #20 hand files and if it is snug at length, the canal is shaped and ready to pack
- If the #20 files is loose at length, proceed to B2. When necessary, continue to B3, gauging after each finishing files with the #25, #30, 40#, #50 respectively.

9. Technique tips

- Using GO-TAPER UNIVERSAL with a torque control, electric endo motor
- Before using GO-TAPER UNIVERSAL, glide path shall be created
- After each rotary files, irrigate, rinsing to clean the canals
- Frequently inspect and clean the flutes of rotary files

10. Cleaning and sterilization

Method of sterilization:

Put this product in a sterilization pack (or foil) and place it on a sterilization tray for autoclave sterilization with reference to the following terms:

1. Sterilize with steam at 134°C (273°F) / 2.1 bars / 18 min.for endodontic instruments and filling.
2. Do not use high pressure steam sterilizer which heats more than 200 degrees Celsius including drying process.
3. When reusing the instrument, wash away foreign substances completely and sterilize.
4. Regarding use of medical cleaning agent, follow the instruction manual by its manufacturer strictly.
5. Dispose the product if damaged or decreased in performance.

11. Warnings

- a. Only skilled dentists are allowed to use.
- b. Be sure to sterilize this product for each use.
- c. Do not use this product except for the dental service and treatment. Use it in accordance with the intended use.
- d. Please confirm that patient has no allergy to Nickel Titanium alloy and other medical reaction before the treatment.
- e. These instruments are intended for use only in dental root canal treatment, users must follow a correct and reasonable way under medical environments.
- f. Please use a rubber dam system to ensure patients' safety when deliver endodontic treatment.
- g. Do not use this product for a patient who indicates sensitization and allergic reaction.

12. Storages and duration of use

- a. Avoid storing at high temperature, humidity and direct sunlight. Keep liquids away. Store it at room temperature.
- b. Do not damage or make a pinhole to packaging materials.
- c. This product is subject to be improved without previous notice. Enforce first-in first-out method for stock management.
- d. Do not store the instrument under a germicidal lamp to avoid deterioration.

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DESIGNED IN FRANCE

Shenzhen SuperLine
Technology Co. Ltd.
Room 1206, Building 4 Chongwen Garden,
Taoyuan Street Nanshan District
Shenzhen 518055 Guangdong China

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Wellkang Ltd. Suite B,