## Ultrasonic Scaler tip

### • Product Model

Please refer to outer packing.

### • Performance and components

This product consists of a head part and a tail part. The tail part is connected to the ultrasonic scaler, and it is replaceable. It is driven by an Ultrasonic Scaler.

### • Scope of application

Used in conjunction with Ultrasonic Scaler, it is suitable for cleaning and shaping of tooth surface, root canal and other parts.

### • Date of manufacture

Please refer to label on packing.

### • Contraindication

- Patients with hemophilia are not allowed to use.
- Patients with cardiac pacemakers are not allowed to use.
- Doctors with cardiac pacemakers are not allowed to use.
- Patients with heart disease, pregnant women and young children should be cautious to use.

### • Operation methods

Recommended power description:

- "Low gear power" means that the device functions normally within 0-30% of the full gear.
- "Medium gear power" means that the device functions normally within 0-60% of the full gear.
- "High gear power" means that the device functions normally within 0-100% of the full gear.

#### 1. *Scaling series:* G mode of Ultrasonic Scaler.

- Recommended service life: 200 operations of all teeth scaling.
- Recommended power: high gear power.

#### 2. *Periodontal treatment series:* used in the Periodontal Treatment mode (P mode) of Ultrasonic Scaler.

- Recommended service life: 50 times.
- Recommended power: medium gear power.

#### 3. *Endo series:* Root Canal Treatment mode (E mode) of Ultrasonic Scaler.

- Recommended service life: 50 root canals.
- Recommended power: low gear power.

### • Cautions

It is normal for the tip to drip and incompletely atomize during operation. The user can adjust to the power level he needs to use.

- ▶ The tip must be tightened to the handpiece and there must be water spray during operation.
- ▶ Tips must be cleaned, disinfected and sterilized before operation.
- ▶ The chemical reagent used in the reprocessing of tip should not be a strong acid or strong alkalin solution with strong corrosiveness.
- ▶ Please replace the tip with a new one when it is broken or worn out.
- ▶ Do not bend or grind the tip.
- ▶ After the golden working tip is exposed to oxidizing solution such as hydrogen peroxide or sodium hypochlorite, please wash it with water in time to avoid damage to the surface of the coating.
- ▶ Do not continue to operate when the work tip is stuck. It should be removed first, otherwise the tip may be broken.

### • Cleaning, disinfection and sterilization

#### Warnings:

The use of strong detergent and disinfectant (alkaline pH>9 or acid pH <5) will reduce the life span of products. And in such cases, the manufacturer takes no responsibility. This device shall not be exposed to high temperature above 138°C.

#### Processing limit:

The products have been designed for a large number of sterilization cycles. The materials used in manufacture were selected accordingly. However, with every renewed preparation for use, thermal and chemical stresses will result in ageing of the products. The maximum number of sterilizations is 300 times.

### 1. Initial processing

#### 1.1 Processing principles

It is only possible to carry out effective sterilization after the completion of effective cleaning and disinfection. Please ensure that, as part of your responsibility for the sterility of products during use, only sufficiently validated equipment and product-specific procedures are used for cleaning/disinfection and sterilization, and that the validated parameters are adhered to during every cycle. Please also observe the applicable legal requirements in your country as well as the hygiene regulations of the hospital or clinic, especially with regard to the additional requirements for the inactivation of prions.

#### 1.2 Post-operative treatment

The post-operative treatment must be carried out immediately, no later than 30 minutes after the completion of the operation. The steps are as follows:

- a. Let the Ultrasonic Scaler works for 20-30 seconds under irrigation mode to flush the handpiece and tip.
- b. Remove the handpiece from the Ultrasonic Periodontal Treatment Device, and rinse away the dirt on the surface of product with pure water (or distilled water/deionized water).
- c. Dry the product with a clean, soft cloth and place it in a clean tray.

#### **Notes :**

The water used here must be pure water, distilled water or deionized water.

- a. Remove the tip from product with endo wrench or torque wrench, and then put the tip and wrench into a clean tray.
- b. Then use soft cloth to dry the product and accessories and put them into a clean tray.

### 2. Cleaning

The cleaning should be performed no later than 24 hours after the operation.

The cleaning can be divided into automated cleaning and manual cleaning. Automated cleaning is preferred if conditions permit.

#### Automated cleaning

- The cleaner is proved to be valid by CE certification in accordance with EN ISO 15883.
- There should be a flushing connector connected to the inner cavity of the product.
- The cleaning procedure is suitable for the product, and the irrigating period is sufficient.

It is recommended to use a washer-disinfector in accordance with EN ISO 15883. For the specific procedure, please refer to the automated disinfection section in the next section "Disinfection".

#### **Notes :**

- a. The cleaning agent does not have to be pure water. It can be distilled water, deionized water or multi-enzyme. But please ensure that the selected cleaning agent is compatible with the product.
- b. In washing stage, the water temperature should not exceed 45°C, otherwise the protein will solidify, and it would be difficult to remove.
- c. After cleaning, the chemical residue should be less than 10mg / L.

### 3. Disinfection

Disinfection must be performed no later than 2 hours after the cleaning phase. Automated disinfection is preferred if conditions permit.

#### 3.1 Automated Disinfection-Washer-disinfector

- The washer-disinfector is proved to be valid by CE certification in accordance with EN ISO 15883.
- Use high temperature disinfection function. The temperature does not exceed 134°C, and the disinfection under the temperature cannot exceed 20 minutes.
- The disinfection cycle is in accordance with the disinfection cycle in EN ISO 15883.

### 3.2 Cleaning and disinfecting steps by using Washer-disinfector

- Carefully place the product into the disinfection basket. Fixation of product is needed only when the product is removable in the device. The products are not allowed to contact each other.
- Use a suitable rinsing adaptor and connect the internal water lines to the rinsing connection of the washer-disinfector.
- Start the program.
- After the program is finished, remove the product from the washer-disinfector, inspect (refer to section "Inspection and Maintenance") and packaging (refer to chapter "Packaging"). Dry the product repeatedly if necessary (refer to section "Drying").

#### Notes :

- Before use, you must carefully read the operating instructions provided by the equipment manufacturer to familiarize yourself with the disinfection process and precautions.
- With this equipment, cleaning, disinfection and drying will be carried out together.
- Cleaning:
  - The cleaning procedure should be suitable for the product to be treated. The flushing period should be sufficient (5-10 minutes). Pre-wash for 3 minutes, wash for another 5 minutes, and rinse it for twice with each rinse lasting for 1 minute.
  - In the washing stage, the water temperature should not exceed 45 °C, otherwise the protein will solidify and it is difficult to remove.
  - The solution used can be pure water, distilled water, deionized water or multi-enzyme solution, etc., and only freshly prepared solutions can be used.
- Disinfection:
  - Direct use after disinfection: temperature  $\geq 90^{\circ}\text{C}$ , time  $\geq 5$  min or  $\text{A0} \geq 3000$ .
  - For the disinfection here, the temperature is  $93^{\circ}\text{C}$ , the time is 2.5 min, and  $\text{A0} > 3000$ .
- Only distilled or deionized water with a small amount of microorganisms ( $<10$  cfu/ml) can be used for all rinsing steps. (For example, pure water that is in accordance with the European Pharmacopoeia or the United States Pharmacopoeia).
- After cleaning, the chemical residue should be less than 10mg / L.
- The air used for drying must be filtered by HEPA.
- Regularly repair and inspect the disinfector.

### 4. Drying

If your cleaning and disinfection process does not have an automatic drying function, dry it after cleaning and disinfection.

#### Methods:

- Spread a clean white paper (white cloth) on the flat table, point the product against the white paper (white cloth), and then dry the product with filtered dry compressed air (maximum pressure 3 bar). Until no liquid is sprayed onto the white paper (white cloth), the product drying is completed.
- It can also be dried directly in a medical drying cabinet (or oven). The recommended drying temperature is  $80^{\circ}\text{C}\sim 120^{\circ}\text{C}$  and the time should be 15~40 minutes.

#### Notes :

- The drying of product must be performed in a clean place.
- The drying temperature should not exceed  $138^{\circ}\text{C}$ .
- The equipment used should be inspected and maintained regularly.

### 5. Inspection and maintenance

In this chapter, we only check the appearance of the product.

- Check the product. If there is still visible stain on the product after cleaning/disinfection, the entire cleaning/disinfection process must be repeated.
- Check the product. If it is obviously damaged, smashed, detached, corroded or bent, it must be scrapped and not allowed to continue to be used.
- If the service time (number of times) of the product reaches the specified service life (number of times), please replace it in time.

### 6. Packaging

Install the disinfected and dried product and quickly package it in a medical sterilization bag (or special holder, sterile box).

#### Notes :

- The package used conforms to ISO 11607.
- It can withstand high temperature of  $138^{\circ}\text{C}$  and has sufficient steam permeability.
- The packaging environment and related tools must be cleaned regularly to ensure cleanliness and prevent the introduction of contaminants.
- Avoid contact with parts of different metals when packaging.

### 7. Sterilization

Use only the following steam sterilization procedures (fractional pre-vacuum procedure\*) for sterilization, and other sterilization procedures are prohibited:

- The steam sterilizer complies with EN13060 or is certified according to EN 285 to comply with EN ISO17665.
- The highest sterilization temperature is  $138^{\circ}\text{C}$ .
- The sterilization time is 18 minutes at  $134^{\circ}\text{C}$  and a pressure of 2.0 bar ~ 2.3 bars.

#### Notes :

- Only products that have been effectively cleaned and disinfected are allowed to be sterilized.
- Before using the sterilizer for sterilization, read the Instruction Manual provided by the equipment manufacturer and follow the instructions.
- Do not use hot air sterilization and radiation sterilization as this may result in damage to the product.
- Please use the recommended sterilization procedures for sterilization. It is not recommended to sterilize with other sterilization procedures such as ethylene oxide, formaldehyde and low temperature plasma sterilization. The manufacturer assumes no responsibility for the procedures that have not been recommended. If you use the sterilization procedures that have not been recommended, please adhere to related effective standards and verify the suitability and effectiveness.  
Fractional pre-vacuum procedure = steam sterilization with repetitive pre-vacuum. The procedure used here is to perform steam sterilization through three pre-vacuums.

### 8. Storage

- Store in a clean, dry, ventilated, non-corrosive atmosphere with a relative humidity of 10% to 93%, an atmospheric pressure of 70KPa to 106KPa, and a temperature of  $-20^{\circ}\text{C}$  to  $+55^{\circ}\text{C}$
- After sterilization, the product should be packaged in a medical sterilization bag or a clean sealing container and stored in a special storage cabinet. The storage time should not exceed 7 days. If it is exceeded, it should be reprocessed before use.

#### Notes :

- The storage environment should be clean and must be disinfected regularly.
- Product storage must be batched and marked and recorded.

### 9. Transportation

- Prevent excessive shock and vibration during transportation, and handle with care.
- It should not be mixed with dangerous goods during transportation.
- Avoid exposure to sun or rain or snow during transportation.

	Date of manufacture		Do not use if package is damaged
	Batch code		Medical Device
	Catalog number		Packaging unit
	Importer		Sterilizable in a steam sterilizer (auto-clave) at the temperature specified
	Non-sterile		Thermo-disinfectable

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**Guilin Woodpecker Medical Instrument Co., Ltd.**
**MedNet EC-Rep GmbH**