GO-TAPER BLUE

1. Specifications

Model : A0, A1, A2, B1, B2, B3, B4, B5
 Length: 19mm, 21mm, 25mm, 31mm
 Packaging: Sterile blister - 6 files

2. Raw materials

NiTi alloy with heat treatment

3. Function

This is the dental instrument with columnar or tapering working part and it possesses circular form in cross cut section and flat end.

4. Recommend Speed and Torque

Model	Length (mm)	Taper	Speed (RPM)	Torque
A0	19	4%	250-350	3.0
A1	21, 25, 31	2%	250-350	3.0
A2	21, 25, 31	4%	250-350	1,5
B1	21, 25, 31	7%	250-350	1,5
B2	21, 25, 31	8%	250-350	3.0
B3	21, 25, 31	9%	250-350	3.0
B4	21, 25, 31	6%	250-350	3.0
B5	21, 25, 31	5%	250-350	3.0

5. Motion

Continuous rotation

6. Precautions for use

- Go-Taper Blue files are supplied sterile and reuse may increase the risk of cross-contamination or breakage. breakage.
- Repeated disinfection and sterilisation may result in an increased risk of file separation and fracture. fracture of the files.
- This instrument should not be immersed in sodium hypochlorite solution.
- Cleaning the instrument: follow the instructions strictly.

7. Selecting the correct GO-TAPER BLUE

- Create straight line access to canal orifice with A0.
- Use #10 and #15 hand files for coronal 2/3. Work with hand files instruments until a smooth glide path is confirmed. Alternatively, glide path files may be used after #10 hand file.
- In the presence of NaOCl, "float" the A1 in the canal and passively "follow" the glide path. Before light resistance is encountered, "brush" and cut dentin to improve straight-line. Always brush away from the furcation.
- Continue shaping with A1 until the depth of the #15 hand file is reached.
- Use the A2, exactly as described for the A1, until the depth of the #15 hand file is reached.
- In the presence of NaOCl, scout the apical 1/3 with #10 and #15 hand files until they are loose at length.
- Establish working length, confirm and verify the presence of glide path in the apical 1/3.
- Use the A1, with a brushing action, until working length is reached.
- Use the A2, with a brushing action, until working length is reached.
- Reconfirm working length, irrigate, recapitulate and re-irrigate, especially in more curved canals.
- Use Finishing File B1, in a "non-brushing" action, with each insertion deeper than the previous insertion until working length is reached. Do not leave the file at working length for longer than one second.
- Gauge the foramen with a #20 hand file. If the instrument is snug at length, the canal is shaped and ready to be obturated.
- If the #20 hand file is loose at length, proceed to the B2 and, when necessary the B3, B4 and B5, with the same non-brushing motion to working length, gauging after each finishing file with #25, #30, #40 or #50 hand files respectively.

8. Sequence

- X-ray film to estimate the working length of canals
- Create straightline access to canal orifice
- In the presence with EDTA, scout the coronal two-thirds with SS hand K-files #10/#15, then gently work K-files until a smooth glide path is confirmed
- Rinsing the canal with NaOCl enough, use A1 and follow the glide path until light resistance, then laterally "brush" and cut dentin on the outstroke

- Shaping with A1 as described until the depth of the #15 hand file is reached
- Use A2, exactly as described for A1, until the depth of the #15 hand file is reached
- In the presence with EDTA, scout the apical one-third with SS hand K-files #10/#15, and gently work them until they are loose at length
- Establish working length, confirm patency and verify a smooth glide path in the apical one-third
- Use A1, as described, until working length is reached
- Use A2, as described, until working length is reached
- Confirm working length
- Use B1 to shape the apical area
- Gauge the shaping canals with #20 hand files and if it is snug at length, the canal is shaped and ready to pack
- If the #20 files is loose at length, proceed to B2. When necessary, continue to B3, gauging after each finishing files with the #25, #30, 40#, #50 respectively.

9. Cleaning and sterilization

Sterilisation method for re-use:

Put this product in a sterilization pack (or foil) and place it on a sterilization tray for autoclave sterilization with reference to the following terms:

1. Sterilize with steam at 134°C (273°F) / 2,1 bars / 18 min.for endodontic instruments and filling.
2. Do not use high pressure steam sterilizer which heats more than 200 degrees Celsius including drying process.
3. When reusing the instrument, wash away foreign substances completely and sterilize.
- 4.Regarding use of medical cleaning agent, follow the instruction manual by its manufacturer strictly.
5. Dispose the product if damaged or decreased in performance.

10. Warnings

- a. Only skilled dentists are allowed to use.
- b. Be sure to sterilize this product for each use.
- c. Do not use this product except for the dental service and treatment. Use it in accordance with the intended use.
- d. Please confirm that patient has no allergy to Nickel Titanium alloy and other medical reaction before the treatment.
- e. These instruments are intended for use only in dental root canal treatment, users must follow a correct and reasonable way under medical environments.
- f. Please use a rubber dam system to ensure patients' safety when deliver endodontic treatment.
- g. Do not use this product for a patient who indicates sensitization and allergic reaction.

11. Storages and duration of use

- a. Avoid storing at high temperature, humidity and direct sunlight. Keep liquids away. Store it at room temperature.
- b. Do not damage or make a pinhole to packaging materials.
- c. This product is subject to be improved without previous notice. Enforce first-in first-out method for stock management.
- d. Do not store the instrument under a germicidal lamp to avoid deterioration.

Assorted package

Do not use if package is damaged

6 pcs/box

Nickel titanium product

Sterilized using irradiation

Stoppers are silicon stoppers

For endo motor use

Sterilize at 134° C

Speed

Continuous rotation

Manufacturer

Authorized representative in the European Community

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Wellkang Ltd.
Shenzhen SuperLine Technology Co. Ltd.