GO-TAPER UNIVERSAL

1. Specifications

Model: A0, A1, A2, B1, B2, B3, B4, B5

Length : 19mm, 21mm, 25mm, 31mm

Package: Sterile blister - 6 files

2. Raw materials

NiTi alloy

3. Function

This is the dental instrument with columnar or tapering working part and it possesses circular form in cross cut section and flat end. GO-TAPER UNIVERSAL files are used aiming at shaping and cleaning the root canal by mechanically powered GO-TAPER UNIVERSAL system.

4. Recommend Speed and Torque

Model	Length (mm)	ISO	Taper	Speed (RPM)	Torque (N.CM)
A0	19	019	4%	250 - 300	5.0
A1	21, 25, 31	018	2%	250 - 300	5.0
A2	21, 25, 31	020	4%	250 - 300	2.0
B1	21, 25, 31	020	7%	250 - 300	2.0
B2	21, 25, 31	025	8%	250 - 300	3.0
B3	21, 25, 31	030	9%	250 - 300	3.0
B4	21, 25, 31	040	6%	250 - 300	3.0
B5	21, 25, 31	050	5%	250 - 300	3.0

5. Motion

Continuous rotation

6. Precautions for use

- Go-Taper Universal files are supplied sterile and reuse may increase the risk of cross-contamination or breakage. breakage.
- Repeated disinfection and sterilisation may result in an increased risk of file separation and fracture. fracture of the files.
- This instrument should not be immersed in sodium hypochlorite solution.
- Cleaning the instrument: follow the instructions strictly.

7. Selecting the correct GO-TAPER UNIVERSAL

- A0- Opener
- A1, A2 work to smooth canal, shaping canal
- B1, B2, B3, B4, B5 are finishing files.
- Normally B2 is enough for most of case.
- If canal is bigger, then use B3, B4, B5.

8. Sequence

- X-ray film to estimate the working length of canals
- Create straightline access to canal orifice
- In the presence with EDTA, scout the coronal two-thirds with SS hand K-files #10/#15, then gently work K-files until a smooth glide path is confirmed
- Rinsing the canal with NaOCl enough, use A1 and follow the glide path until light resistance, then laterally “brush” and cut dentin on the outstroke
- Shaping with A1 as described until the depth of the #15 hand file is reached
- Use A2, exactly as described for A1, until the depth of the #15 hand file is reached
- In the presence with EDTA, scout the apical one-third with SS hand K-files #10/#15, and gently work them until they are loose at length
- Establish working length, confirm patency and verify a smooth glide path in the apical one-third
- Use A1, as described, until working length is reached
- Use A2, as described, until working length is reached
- Confirm working length
- Use B1 to shape the apical area
- Gauge the shaping canals with #20 hand files and if it is snug at length, the canal is shaped and ready to pack
- If the #20 files is loose at length, proceed to B2. When necessary, continue to B3, gauging after each finishing files with

the #25, #30, 40#, #50 respectively.

9. Technique tips

- Using GO-TAPER UNIVERSAL with a torque control, electric endo motor
- Before using GO-TAPER UNIVERSAL, glide path shall be created
- After each rotary files, irrigate, rinsing to clean the canals
- Frequently inspect and clean the flutes of rotary files

10. Cleaning and sterilization

Sterilisation method for re-use:

Put this product in a sterilization pack (or foil) and place it on a sterilization tray for autoclave sterilization with reference to the following terms:

1. Sterilize with steam at 134°C (273°F) / 2,1 bars / 18 min.for endodontic instruments and filling.
2. Do not use high pressure steam sterilizer which heats more than 200 degrees Celsius including drying process.
3. When reusing the instrument, wash away foreign substances completely and sterilize.
4. Regarding use of cleaning agent, follow the instruction manual by its manufacturer strictly.
5. Dispose the product if damaged or decreased in performance.

11. Warnings

- a. Only skilled dentists are allowed to use.
- b. Be sure to sterilize this product for each use.
- c. Do not use this product except for the dental service and treatment. Use it in accordance with the intended use.
- d. Please confirm that patient has no allergy to Nickle Titanium alloy and other medical reaction before the treatment.
- e. These instruments are intended for use only in dental root canal treatment, users must follow a correct and reasonable way under medical environments.
- f. Please use a rubber dam system to ensure patients' safety when deliver endodontic treatment.
- g. Do not use this product for a patient who indicates sensitization and allergic reaction.

12. Storages and duration of use

- a. Avoid storing at high temperature, humidity and direct sunlight. Keep liquids away. Store it at room temperature.
- b. Do not damage or make a pinhole to packaging materials.
- c. This product is subject to be improved without previous notice. Enforce first-in first-out method for stock management.
- d. Do not store the instrument under a germicidal lamp to avoid deterioration.

Assorted package

Do not use if package is damaged

6 pcs/box

Nickel titanium product

Sterilized using irradiation

Stoppers are silicon stoppers

For endo motor use

Sterilize at 134° C

Speed
250 -
350 rpm

Continuous rotation

Manufacturer

Authorized representative in the
European Community

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Wellkang Ltd.
Shenzhen SuperLine Technology Co. Ltd.