1. Specifications

Total length: 21, 25 and 31 mm

Working part length: 16 and 21 mm

Package: Sterile blister - 6 files

3 assorted kits available: ISO N°10-20 for working length 16 mm and 21 mm

ISO N°25-20 for working length 16 mm

2. Raw materials

NiTi alloy

3. Function

Endodontic sequence of constant taper, in Nickel Titanium (NiTi), for the shaping and cleaning of the root canal and the widening canals cleaning.

For standard and difficult cases.

Fits to all endodontic motors and contra-angles with rotary movement.

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Continuous rotation

6. Precautions for use

- Safety equipment (gloves, glasses, rubber dam) is required for patients and self-protection.
- These instruments can only be used after sterilization, please using these filed under medical environment.
- Please use GO-TWO system in an appropriate case and in a correct way. Remove debris out frequently from screwing-in problems
- Normally shaping the canal with EDTA for entering the canal smoothly and irrigate sufficiently after removing debris created by shaping the root canal.
- If you come across a situation that it's difficult to enter the deeper part of the canal, please stop and check both morphology of the canal and whether there are signs of deformation of files or not, then confirm the orifice is expanded wide enough for subsequent use of shaping files.
- GO-TWO system files are supplied sterile and reuse may increase the risk of cross-contamination or breakage.

7. Selecting the correct GO-TWO

- Basic sequence: 10/04, 15/05, 20/06, 25/06
- Shaping for a wider canal: 30/05, 35/04, 40/04, 45/04, 50/04, 60/04
- Shaping for a hot filling: 25/07, 30/06, 35/06, 40/06

8. Sequence

- After opening the chamber and locating all canals, fill the access cavity with sodium hypochlorite.
- Take a K10 file and insert it passively into the canal, rotating it 1/8th of a turn to the left and right. Move forward without applying pressure and stop when the file stops moving. Note this length.
- Take your first instrument from the series (GO-TWO with purple ring) and prepare the channel to this length. To do this, insert the instrument in motion into the canal and move forward passively towards the apex. Move it up brushing the canal wall. Remove the file and clean the flute. Repeat this action until you have reached the desired length. Then do the same with the second instrument of the series (white ring). You have now prepared the coronal 2/3 of the canal.
- Irrigate abundantly.

- Take again your K10 file and move forward passively like before. In the majority of cases, the file will reach the working length passively. Measure this length using an apex locator. If at this point you cannot reach the working length, do not force it. Remove the K10 file and prepare the canal portion obtained using your purple and white GO-TWO instruments, then repeat this step.
- Prepare the canal to the working length using the whole series. Start with the purple GO-TWO, then use the white, yellow, and finally the red instrument. Irrigate well between each instrument. The mechanical preparation of the tooth is thus completed.

9. Technique tips

- X-ray is necessary before retreatment.
- Analysis the root canal carefully.
- Obtain straight line access to the canal.
- Heat will be helpful for filling materials removal.
- Go-Two files are sterilizable
- Smooth channel
- The basic sequence: 10/04 is the first file at the working length with a «brushing» movement. As soon as the working length is reached, move on to the next instruments in the basic sequence: 15/05, 20/06, 25/06. All files are brought to the working length.
- Preparation with a larger file
- If the hand file No. 30 can be brought to 1 mm less than the working length, the apical area should not need to be widened further.
- However, if the No. 30 hand file can be inserted up to the working length, additional preparation is necessary. Please use the following sizes: 30/05, 35/04, 40/04, 45/04, 50/04, 60/04.

10. Cleaning and sterilization

Sterilisation method for re-use:

Put this product in a sterilization pack (or foil) and place it on a sterilization tray for autoclave sterilization with reference to the following terms:

1. Sterilize with steam at 134°C (273°F) / 2,1 bars / 18 min for endodontic instruments and filling.
2. Do not use high pressure steam sterilizer which heats more than 200 degrees Celsius including drying process.
3. When reusing the instrument, wash away foreign substances completely and sterilize.
4. Regarding use of cleaning agent, follow the instruction manual by its manufacturer strictly.
5. Dispose the product if damaged or decreased in performance.

11. Warnings

- a. Only skilled dentists are allowed to use.
- b. Be sure to sterilize this product for each use.
- c. Do not use this product except for the dental service and treatment. Use it in accordance with the intended use.
- d. Please confirm that patient has no allergy to Nickel Titanium alloy and other medical reaction before the treatment.
- e. These instruments are intended for use only in dental root canal treatment, users must follow a correct and reasonable way under medical environments.
- f. Please use a rubber dam system to ensure patients' safety when deliver endodontic treatment.
- g. Do not use this product for a patient who indicates sensitization and allergic reaction.

12. Storages and duration of use

- a. Avoid storing at high temperature, humidity and direct sunlight. Keep liquids away. Store it at room temperature.
- b. Do not damage or make a pinhole to packaging materials.
- c. This product is subject to be improved without previous notice. Enforce first-in first-out method for stock management.
- d. Do not store the instrument under a germicidal lamp to avoid deterioration.

6 pcs/box

Continuous rotation

Sterilized using irradiation

Nickel titanium product

For endo motor use

Stoppers are silicon stoppers

250 rpm

Speed

Sterilize at 134° C

Do not use if package is damaged

Manufacturer

Importer

Authorized representative in the European Community

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Wellkang Ltd.
Shenzhen SuperLine Technology Co. Ltd.