access

EASY FILL

- (en) Obturation system for root canal filling
- (fr) Système d'obturation du canal radiculaire
- (it) Sistema di otturazione del canale radicolare
- (es) Sistema de obturación del conducto radicular
- (pt) Sistema de obturação de canais radiculares
- (de) Obturationssystem für Wurzelkanalfüllungen
- (nl) Obturatiesysteem voor wortelkanaalvulling

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EN

1. Product introduction

- 1.1 Structural composition
- 1.2 Package contents
- 1.3 Scope of application
- 1.4 Contraindications
- 1.5 Equipment safety classification
- 1.6 Main technical parameters and environment of use

2. Product installation and function description

- 2.1 Schematic diagram of the components
- 2.2 Power adapter installation
- 2.3 Assembly and disassembly of the gutta percha injection needle
- 2.4 Use of each function button
- 2.5 Functional operation description
- 2.6 Loading the gutta percha bar

3. Instructions for use of accessories

- 3.1 Selection of the gutta percha injection needle
- 3.2 Using the silver needle bending device
- 3.3 Gutta percha bar selection

4. Precautions

5. Troubleshooting

6. Cleaning, disinfection and sterilization

- 6.1 Cleaning of the charging base and the central unit
- 6.2 Cleaning and disinfection of the gutta percha injection needle, silver needle bending device and heat insulation protection cover
- 6.3 High-temperature and high-pressure sterilization

7. Storage, maintenance and transportation

- 7.1 Storage and maintenance
- 7.2 Transportation

8. Environmental protection

9. After-sales service

10. Electromagnetic compatibility statement

11. Symbol description

1. Product introduction

The **Easy Fill** obturation system is suitable for root canal filling during root canal treatment. It is used to heat and soften the gutta-percha, so that the gutta-percha can be fully filled into the root canal.

This equipment features:

- Easy and convenient electric pushing of the gutta percha.
- Wireless «gun-style» design for comfortable grip.
- Clear, simple and easy-to-operate LED screen.
- Accurate temperature control, with four preset temperatures that can be selected: 100°C, 150°C, 180°C, 200°C.
- Safety protection mechanism, which will automatically shut the unit down after three minutes without operation.

1.1 Structural composition

The device is composed of a central unit, a charging base, a power adapter, a gutta percha injection needle and a heat insulation protection cover; main parts of the central unit are the display screen, the lithium battery and the handheld part.

1.2 Package contents

Name	Quantity	Name	Quantity
Central unit	1	Silver needle bending device	1
Charging base	1	Heat insulation protection cover	3
Power adapter	1	Gutta percha injection needle	6
Dual-headed Type-C charging cable	1	Instruction	1
Type-C charging cable	1		

1.3 Scope of application

Suitable for root canal filling in root canal therapy.

1.4 Contraindications

- Prohibited for patients with known allergy to natural latex, stainless steel, silver, copper and other metals.
- Prohibited for patients with hemophilia.
- Prohibited for patients with pacemakers.
- Prohibited for use by physicians with pacemakers.
- Use with caution in patients with heart disease, pregnant women and young children.

1.5 Equipment safety classification

- Classification by operation mode: short-time operation.
- Classification by electric shock protection type: class II equipment.
- Classification by electric shock protection degree: type-B applied part.
- Protection against ingress of fluids: ordinary equipment (IPX0).
- Not to be used with a flammable anesthetic gas mixed with air or a flammable anesthetic gas mixed with oxygen or nitrous oxide.

1.6 Main technical parameters and environment of use

Size	Central unit: 150 x 124 x 30 mm
Weight	About 186 g
Power supply mode	Lithium battery, DC3.7V, 2600 mAh
Power adapter	Input: AC100-240V, 50/60Hz 0.4A Output: DC5V/2A
Fluid penetration protection	IPX0
Classification by electric shock protection type	Class II, internal power supply equipment
Classification by electric shock protection degree	Class B
Optional temperature	100°C, 150°C, 180°C, 200°C

Power dissipation	<3W
Service environment	Temperature requirements: 5°C ~ 40°C Humidity requirements: 30% ~ 75% Atmospheric pressure: 70 kPa ~ 106 kPa
Storage / transport environment:	Temperature requirements: -20°C ~ 55°C Humidity requirements: 10% ~ 93% Atmospheric pressure: 500 hPa ~ 1060 hPa

2. Product installation and function description

2.1 Schematic diagram of components

1. Power button
2. LED display screen
3. Setting button
4. Central unit
5. Charging base
6. Dual-headed Type-C charging cable

7. Type-C charging cable
8. Silver needle bending device
9. Power adapter
10. Gutta percha injection needle
11. Heat insulation protection cover

2.2 Power adapter installation

This product uses contact charging. When the power is insufficient, please place the unit on the charging base in time for charging.

The power adapter and the base should be connected to both ends of the Type-C charging cable as shown in Figure (a) (b), and the dual-headed Type-C charging cable can be connected as shown in Figure (c) for an output terminal.

When the central unit is correctly connected to the charging base, the charging base display light changes from blue to flashing green; after full charge, the light turns to steady green.

(a)

(b)


(c)

2.3 Assembly and disassembly of the gutta percha injection needle


2.3.1 Installing the gutta percha injection needle

Please remove the protective plastic head before the first use. As shown in the figure, select the required gutta percha injection needle and attach it by turning it clockwise.

2.3.2 Installing the heat insulation protection cover

Heat insulation protection cover on the gutta percha injection needle

As shown in the figure, please place the heat insulation protection cover on the gutta percha injection needle to prevent burns from working heating.

2.3.3 Removing the gutta percha injection needle

After use, when the device is cooled to room temperature, carefully remove the heat insulation protection cover, then rotate the gutta percha injection needle counterclockwise to remove it, as shown in the figure, and put it into the gutta percha injection needle holder after disinfection.

2.4 Use of each function button

2.4.1 Power button

When the device is off, long press the power button to switch it on.
When the device is on, long press the power button to switch it off.

2.4.2 Setting button

With the device switched on, short press the setting button to select the four preset temperatures, as shown in Figure (a).

Long press the setting button and then short press it to switch between the three options shown in Figure (b).

(a)

(b)


2.4.3 Trigger

After determining the temperature required in the main interface, press the trigger short for heating. After heating, press the trigger long to push the gutta percha. When changing mode, press the trigger to confirm the current option and execute.

2.5 Functional operation description

2.5.1 Main interface

- After entering the main interface, the default temperature value (T), the remaining gutta percha amount (GP) and the actual power of the battery are displayed on the display screen as shown in the figure.

- Battery icon:

The figure below shows the battery power display icon, which is dynamic during battery charge.

				
80-100% electric charge	60-80% electric charge	40-60% electric charge	20-40% electric charge	: 0-20% charge, the battery is depleted. Please charge the device immediately.

Caution : This icon indicates the remaining battery power. During gutta percha filling, the battery power will decrease.

2.5.2 Operation Interface



Long press the setting button to enter the operation interface: ① RELOAD GP (automatic residual gutta percha removal and reset), ② EXIT return to the main interface). Press the setting button to select EXIT (return

to the main interface) and press the trigger to return to the main interface. ③ **RESET** (quick reset). Press the setting button to select RESET (quick reset). Press the trigger to start the quick reset program. After the quick reset, the display screen displays RESET COMPLETE.

① **RELOAD GP** (automatic residual gutta percha removal and reset). Short press the setting button to select RELOAD GP (automatically remove the residual gutta percha and reset). Press the trigger to start the program to automatically remove the remaining gutta percha and reset. After the display screen shows RELOAD GP, the dental gutta percha bar can be installed as described in section 2.6.

② **EXIT** (return to the main interface). Press the setting button to select EXIT (return to the main interface) and press the trigger to return to the main interface.


③ **RESET** (quick reset). Press the setting button to select RESET (quick reset). Press the trigger to start the quick reset program. After the quick reset, the display screen displays RESET COMPLETE.

• 2.5.3 Heating interface

Short press the trigger for heating, if you need to stop heating, short press the power button.

The contents displayed on the heating interface as shown in the Figure (a) are: heating state (HT), real-time temperature and remaining gutta percha quantity progress bar. The indicator flashes when heating starts, and it is always on after heating. After heating, the contents displayed on the heating interface as shown in Figure (b) are: heating completion state (OK), preset temperature and remaining gutta percha quantity progress bar.

(a)

(b)

• 2.5.4 Push interface

After heating, long press the trigger to push the gutta percha. The gutta percha push interface, shown in the figure, includes the remaining gutta percha quantity progress bar, the heating completion state (OK), the preset temperature and the direction indication arrow. When the gutta percha bar is exhausted, reset starts automatically.

2.6 Loading the gutta percha bar

Before loading, you need to confirm that the device is in the state of «reset complete» (in the case of newly purchased products, the gutta percha bar can be directly loaded without using the «reset» state). The gutta percha bar must first be inserted into the gutta percha injection needle, then attach the gutta percha injection needle to the device as described in section 2.3.

Attention :

Before loading the gutta percha bar, please confirm that the gutta percha remaining in the heating cavity has been removed. If not, please perform the RELOAD GP (automatic removal and reset) procedure on the operation interface (refer to 2.5.2 for details) before loading. If the device is automatically reset, the gutta percha can be loaded.

3. Instructions for use of accessories


3.1 Selection of the gutta percha injection needle

Before use, the gutta percha injection needle should be installed on the central unit of the gutta percha filling device. Different sizes of gutta percha injection needle are provided, as shown in the following table:

Model	Size	Length
S23G-24	23G	24mm
S25G-24	25G	24mm
S25G-28	25G	28mm

3.2 Using the silver needle bending device

One end of the device shown in Figure (a) can be used to bend the needle with a radius of 5mm or 3mm (Figure b), and the other end can be used to rotate it 360° (Figure c).

a)

b)

c)

3.3 Gutta percha bar selection

Gutta percha bar size: diameter 2.5-2.8 mm, length less than 14 mm.

4. Precautions

- As it is the case with all electronic devices, this instrument causes electromagnetic interference and should not be used in patients with pacemakers.
- During operation and cleaning, pay attention to the direction and position of the gutta percha injection needle to avoid injuring patients and operators.
- Pay full attention to whether the gutta percha injection needle is loose, hot or makes noise. If you detect any abnormal phenomenon, please stop using the device and contact the local dealer or manufacturer.

- Keep the device safe from impacts and especially from dropping.
- When removing the gutta percha injection needle and heat insulation protection cover, turn off the power first in order to avoid any accidental starting and injury that might result from touching the trigger accidentally.
- When the battery power icon shows low power, please charge in time using the original power adapter.
- Please use the original accessories, those from other manufacturers may not work with this equipment or cause equipment damage.
- Please do not disassemble or repair the equipment without permission, otherwise the warranty will no longer be valid.
- This equipment should not be placed in wet places or be in contact with any liquid.
- Do not expose the equipment to a direct or indirect heat source. The equipment should be operated and kept in a safe environment.

5. Troubleshooting

Fault	Possible cause	Solution
The device cannot be switched on	1. Insufficient battery power. 2. Short circuit of the charging interface. 3. The hot melt gutta percha filling device is damaged.	1. Charge the device connecting it to a power supply. 2. Remove the foreign bodies at the interface, wipe it clean and reconnect it to the power supply. 3. Contact the local dealer or our company.
The gutta percha does not flow out from the injection needle	1. The gutta percha quantity progress bar has reached the end, indicating the gutta percha bar has been used up. 2. The gutta percha injection needle is damaged or clogged.	1. Turn on the gutta percha load function manually. 2. Replace the gutta percha injection needle.
The gutta percha injection needle cannot be removed	The melted gutta percha material has cooled down in the cavity and hardened	1. Turn on the power switch, press the trigger, and pull the needle out when the heating indicator is on. 2. Contact the customer service.
With the power adapter plugged in, the device is not charged	1. Not correctly plugged in. 2. Power supply is damaged or the specifications do not match. 3. The charging base contacts are dirty.	1. Pull out and reconnect. 2. Contact your dealer. 3. Wipe the charging base contacts with alcohol.
With the battery fully charged, the service cycle becomes shorter	Reduced battery capacity.	1. Send the device to the maintenance department for maintenance. 2. Contact the customer service.

6. Cleaning, disinfection and sterilization

6.1 Cleaning of the charging base and the central unit

The charging base and the surface of the central unit can be wiped with a soft cloth dipped in a small amount of neutral detergent or alcohol.

6.2 Cleaning and disinfection of the gutta percha injection needle, silver needle bending device and heat insulation protection cover

Step	Process
1. Rinse	Wash the gutta percha needle, silver needle bending device, heat insulation protection cover with running water for 2 min to remove the surface contaminants.
2. Wipe	Wet a soft clean cloth with a cleaning agent and thoroughly wipe the surface of the product 5 times. Replace the soft cloth after each wipe. If there are still visible pollutants remaining, wipe repeatedly until there are no visible pollutants.
3. Brush wash	Brush the product thoroughly with an instrument brush for 3 min.
4. Soak	Immerse the product in the cleaning agent for 5 min.

5. Rinse	Wash the product with purified water for 2 min to remove any residual detergent from its surface.
6. Dry	Use a dry soft absorbent cloth to eliminate any residual water trace on the product surface.

6.3 High-temperature and high-pressure sterilization

Sterilization	Put the gutta percha injection needle, silver needle bending device and heat insulation protection cover into a disposable sterilization bag and carry out sterilization at a temperature of 134°C during 18 minutes.
---------------	---

7. Storage, maintenance and transportation

7.1 Storage and maintenance

This equipment should be kept away from any source of heat and installed or kept in a cool, dry and ventilated place.

The equipment should not be stored with toxic, corrosive, flammable and explosive articles.

The equipment should be stored in an environment with a relative humidity of 10% -93%, an atmospheric pressure of 500 hPa ~ 1060 hPa and a temperature of -20°C ~ + 55°C.

When the equipment is not in use, turn off the power switch and unplug the charging base.

Using the equipment will shorten the battery life. Please charge in time when the power is low.

When the equipment is not in use for a long time, it should be charged once a month for 1 hour.

7.2 Transportation

Prevent excessive impact and vibration during transportation.

During transportation, the product should not be mixed with dangerous goods.

Avoid exposure to the sun or rain or snow during transportation.

8. Environmental protection

The equipment does not contain any harmful ingredients and can be treated or disposed of in accordance with relevant local policies.

9. After-sales service

As of the date of sale, if the equipment does not work normally due to quality problems, our company shall be responsible for the maintenance with the use of the warranty card. The warranty period and warranty scope shall refer to the product warranty card. This product does not contain parts suitable for self-repair, any maintenance should be carried out by designated professionals or a special maintenance shop. If you need to repair the equipment components, the circuit diagram, component list, correction details and other information can be obtained by contacting the manufacturer.

10. Electromagnetic compatibility statement

The Easy Fill obturation system meets the requirements of the YY0505 standard on electromagnetic compatibility and must be installed and used according to the EMC information specified in this manual. Portable and mobile RF communication equipment may affect this device.

Except for cables (transducers) sold as spare parts for internal components, the use of accessories and cables (transducers) other than those specified for the device may result in increased emissions or decreased immunity of the equipment or system.

The equipment or system should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the device should be observed to verify normal operation in the configuration in which it will be used.

10.1 Guidance and manufacturer's declaration - electromagnetic emissions

Guidance and manufacturer's declaration - electromagnetic emissions		
The Easy Fill obturation system is intended for use in the electromagnetic environment specified below. The customer or the user of the device should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - Guidance
Radio-frequency emissions GB-4824	Group 1	The Easy Fill obturation system uses RF energy only for its internal function. Therefore, its RF emissions are low and are not likely to cause any interference in nearby electronic equipment.

Émissions de radiofréquences GB-4824	Class B	The Easy Fill obturation system is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage network that supplies buildings used for domestic purposes
Radio-frequency emissions GB-4824	Class A	
Harmonic emissions GB-17625.1	Complies	

10.2 Guidance and manufacturer's declaration - electromagnetic immunity

Guidance and manufacturer's declaration - electromagnetic immunity

The Easy Fill obturation system is intended for use in the electromagnetic environment specified below. The customer or the user of the device should ensure that it is used in such an environment.

Immunity test	IEC-60601 test level	Compliance level	Electromagnetic environment - Guidance
Electrostatic discharge GB/T 17626.2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst GB/T 17626.4	± 2kV for power supply lines	± 2kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge GB/T 17626.5	± 1 kV line to line ± 2 kV line to earth	± 1 kV line to line	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines GB/T 17626.11	<5 % UT (> 95% dip in UT) for 0.5 cycle 40 % UT (60% dip in UT) for 5 cycles 70 % UT (30% dip in UT) for 25 cycles <5 % UT (> 95% dip in UT) for 5 sec.	<5 % UT (> 95% dip in UT) for 0.5 cycle 40 % UT (60% dip in UT) for 5 cycles 70 % UT (30% dip in UT) for 25 cycles <5 % UT (> 95% dip in UT) for 5 sec.	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Easy Fill obturation system requires continued operation during power mains interruption, it is recommended that the Easy Fill obturation system be powered from an uninterruptible power supply or battery.
Power frequency (50/60Hz) magnetic field GB/T 17626.8	3 A/m	3 A/m (50/60Hz)	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: UT is the AC mains voltage prior to application of the test level.

10.3 Guidance and manufacturer's declaration - electromagnetic immunity

Guidance and manufacturer's declaration - electromagnetic immunity

The Easy Fill obturation system is intended for use in the electromagnetic environment specified below. The customer or the user of the device should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - Guidance
---------------	----------------------	------------------	--

Conducted RF GB/T 17626.6	3 V (valid value) 150 kHz~80MHz	3 V (valid value)	Portable and mobile RF communications equipment should be used no closer to any part of the Easy Fill obturation system including cables than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2$ $d = 1.2\sqrt{P}$ 80MHz~800MHz $d = 2.3\sqrt{P}$ 800MHz~2.5GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a) should be less than the compliance level in each frequency range, b) interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF GB/T 17626.3	3 V/m 80MHz~2.5GHz	3 V/m	

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

Field strengths from fixed transmitters, such as base stations for radio (cellular/ cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Easy Fill obturation system is used exceeds the applicable RF compliance level above, the Easy Fill obturation system should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Easy Fill obturation system. Over the frequency range 150 kHz to 80MHz, field strengths should be less than 3 V/m.

10.4 Recommended separation distances between portable and mobile RF communications equipment and Easy Fill

The Easy Fill obturation system is intended for use in an environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

	Separation distance according to frequency of transmitter in meter		
Rated maximum output power of transmitter (W)	150 kHz~80MHz $d = 1.2\sqrt{P}$	80MHz~800MHz $d = 1.2\sqrt{P}$	800MHz~2.5GHz $d = 2.3\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80MHz and 800MHz, the separation distance for the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

11. Symbol description

	Consult instruction for use		Caution
	Type B applied part		Medical device
	The product complies with the WEEE directive, and the equipment must be treated as municipal solid waste when disposed of		Product serial number
	Manufacturer		Date of manufacture
	Indoor use		Upward
	Keep dry		Fragile, handle with care
	UDI		

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