GUTTA-PERCHA CONDENSER

1. Description of the product

Our GUTTA-PERCHA CONDENSERS are dental canal roots instruments designed to mechanically condense Gutta Percha in the root canal, once preparation, cleaning and disinfection are finished. It is used with a low-speed contra-angle. These are non-sterile instruments.


2. Product Structure and materials

Our GUTTA-PERCHA CONDENSERS are composed of an active part and a shank. There are made of stainless steel. The rotation of this instrument against the Gutta-Percha plasticizes the material by heat produced from the friction. In this plastic state the Gutta-Percha is advanced towards the apex and condensed.

It is a standard instrument with 2% taper, which is type 1 classification according to ISO 3630-1:2019.

Main materials:

- a. Active part: Stainless Steel
- b. Shank: Brass C3604

3. Before Using

WARNINGS

- Only skilled dentists are allowed to use.
- Ensure to sterilize this product before each use.
- Only use this product for dental care and treatment. Use it in accordance with its intended purpose.
- The GUTTA-PERCHA CONDENSERS must not be used in curved canals, as in this clinical situation there is a high risk of instrument breakage.
- The use of GUTTA-PERCHA CONDENSERS is not risk-free. The main danger consists in a possible penetration of the Gutta-Percha into the periapical region. The following reasons may be responsible for such a case: the use of a too-small or too-fine master cone, a too-small GUTTA-PERCHA CONDENSERS, excessive speed or axial pressure, or too-deep insertion of the GUTTA-PERCHA CONDENSERS in the root canal. Even a small amount of Gutta-Percha in the periapical region may cause severe clinical reactions. There is no resorption expected from the periapical macrophages.
- For multicuspid mandibular teeth, overfilling reaching the inferior dental nerve may provoke irreversible damage, such as dysaesthesia, paraesthesia, anaesthesia, or paralysis with a possible labiodynisis.
- Never use the Gutta-Percha Condenser in counterclockwise rotation (risk of it getting stuck in the canal or in the bone).
- Do not insert the tip of the instrument into the canal to a distance of less than 2 mm from the apex.
- Do not apply strong resistance when the Gutta-Percha Condenser is backing out of the canal.

Contraindications and prohibitions

1. Do not use this product on a patient who indicates sensitization or an allergic reaction.
2. The Gutta-Percha Condenser is not recommended for use in patients who have a known allergy to nickel. Use of this product in these patients may cause difficulty breathing, swollen face or eyes, hives or rash.
3. Should any of these symptoms occur, discontinue use immediately and seek medical advice.

4. Re-Processing (based on ISO 17664: 2017)

Our Dental Root Canal Instruments are single-use devices and are not approved for re-use.

5. Specifications and Method of Use

- The Specifications of ACCESS GUTTA-PERCHA CONDENSERS are shown in the table below (based on ISO 3630-1&ISO3630-5):

Size	L ₃ (mm)	Tolerance	L _n (mm)	Tolerance	L _t (mm)	Tolerance	ΦD (mm)	Φd ₃ (mm)	Tolerance	Φd _n (mm)	Tolerance	Color
#25	3	±0.50	17	±0.50	L21 L25 L31	±0.50	0.25	0.31	±0.02	0.8	±0.02	Red
#30							0.30	0.36				Blue
#35							0.35	0.41				Green
#40							0.40	0.46				Black
#45							0.45	0.51		1.00		White
#50							0.50	0.56				Yellow
#55							0.55	0.61				Red
#60							0.60	0.66		1.20		Blue
#70							0.70	0.76				Green
#80							0.80	0.86				Black

- In order to obtain sufficient heat to plastify the Gutta-Percha, the GUTTA-PERCHA CONDENSERS must be used at a low speed with a high torque instrument capable of rotating at 8'000 rpm.
- For a complete filling of the root canal system, bring the instrument to a distance of 2 mm from the apex of the canal which has already been prepared. The penetration depth of the GUTTA-PERCHA CONDENSERS is checked by using the graduation marked on the shaft of the instrument or by using a silicone-stop.
 - Establish the working length and apical stop of the root canal.
 - Insert the Gutta-Percha point and the GUTTA-PERCHA CONDENSERS into the canal and press the GUTTA-PERCHA CONDENSERS against the Gutta-Percha.
 - Descend the GUTTA-PERCHA CONDENSERS into the canal, rotating in a clockwise direction at 8'000 rpm. Do not penetrate the GUTTA-PERCHA CONDENSERS to a distance of less than 2 mm from the apex.
 - Withdraw the GUTTA-PERCHA CONDENSERS out of the canal, without applying strong resistance.

6. Precautions for use

- To prevent infection, clean and disinfect the product (See section 7) and ensure sterilization is completed before using.
- Choose the most appropriate type for each case and follow the general method.
- Before use, make sure that the instruments outside the oral cavity show no deformations, scratches and cracks.
- If the head of product is thin, long or large, there is a risk of breakage or twisting. Therefore, be sure to avoid using unreasonable angle and excessive pressure.
- Use dental dam or similar protective devices to avoid accidental ingestion and dropping.
- Do not use this instrument for any purpose other than the applications listed above.
- Only for use by dentists.
- This product should be treated as medical waste when disposed.
- Dispose of the product if it is damaged or contaminated.
- After using it, wash it with medical cleaning agent and brush, then remove foreign substances such as adherent body fluids and tissue residues.
- Place the product on a stand when cleaning with an ultrasonic cleaner to avoid deterioration of the cutting part.
- Use this product with great care to avoid puncturing fingers due to sharp edges.
- This product may be subject to corrosion if immersed for a long time in NaOCl, EDTA, or similar solutions.

7. Cleaning and Disinfection (based on ISO 17664)

General Recommendation

- Use cleaning agents suitable for medical instruments, such as multi-enzyme cleaning agents. When using, strictly follow the manufacturer's instructions regarding concentration, temperature, contact time, and shelf life of the solution. Do not use hydrogen peroxide (H₂O₂) solution.
- For your own safety, wear personal protective equipment (gloves, glasses, masks).
- The user is responsible for the sterilization or disinfection of the product during the first cycle, as well as for the use of damaged or insufficiently cleaned equipment after sterilization.
- In the final rinsing step, deionized water must be used regardless of whether an automatic washer sterilizer or manual cleaning method is used. Tap water can also be used for other rinsing steps.
- Only immerse patient-contacting operative part of the instrument in NaOCl solution with a concentration not exceeding 5%.
- Before or during disinfection or cleaning, avoid equipment to dry out. Dried biological materials may be difficult to remove.

- Only use the appropriate support of the device for post-processing.
- Do not use labeling systems or identification marks directly on the device.
- Only use washer-disinfectors approved in accordance with EN ISO 15883 and regularly maintained and validated.

Operation Steps

Disassembling: Remove the stopper from the instrument and dispose of it

Rinsing: Rinse extensively (at least 1 minute) under flowing deionized water (ambient temperature). While rinsing, use a soft brush (made of nylon, polypropylene, acrylic) for pre-cleaning until visible impurities are removed.

Cleaning and Disinfection

Manual Cleaning assisted by an ultrasonic device:

- Manual cleaning with ultrasonic equipment
- Place the instrument in a kit, stand or container (made of stainless steel, polypropylene or titanium).
- Immerse in a detergent solution with cleaning properties, if appropriate, soak for at least 15 minutes with ultrasonic assistance.
- Flushing: Perform an abundant flushing (at least 1 minute) under flowing deionized water temperature 20°C ~ 40°C.
- Drying: Dry with a disposable non-woven fabric or hot air dryer not exceeding 110°C.

Automated Cleaning using a cleaning and disinfecting device

- Place the instrument in a kit, stand or container (made of stainless steel or titanium).
 - Perform the defined cycle with detergent solution for at least 5 minutes in the washer-disinfector with temperature 20°C ~ 40°C.
-

Disinfection using a washer-disinfector device

- Place the instrument in a kit, stand or container (made of stainless steel or titanium).
- Perform the defined cycle with mild neutral enzyme cleaning agent solution for at least 5 minutes in the washer-disinfector with temperature >90°C, A0 >3000.

Note:

1. Discard any instruments with obvious defects (broken, bent, etc.).
 2. When the instruments are placed in the cleaning kit, support or container, avoid any contact between them.
 3. Follow the instructions and concentration provided by the detergent solution manufacturer (see also general recommendations).
 4. Follow the instructions of the washer-disinfector and verify the success criteria after each cycle is reached according to the manufacturer's instructions.
 5. The final rinse step should use deionized water. For other steps, follow the water quality defined by the manufacturer. Place the devices in a kit, support, or container (made from stainless steel or titanium) to avoid any contact between devices or posts.
-

Rinsing

- Abundant rinsing (at least 1 min) under running water (ambient temperature).
1. Use deionized water for rinsing.
 2. If the previously used cleaning solution contains a corrosion inhibitor, it is recommended to do the rinsing step just before starting the autoclaving.
-

Drying:

Devices should be thoroughly dried before inspection and packaging.

1. Dry on a single-use non-woven cloth, or with a hot air dryer at a temperature not exceeding 110°C.
 2. Devices should be dried until visual traces of moisture are eliminated.
 3. Pay particular attention to effectively dry joints or cavities within a device.
-

Inspection

- Inspect the functionality of the devices.
 - Inspect devices and sort out any defective ones.
1. Dirty devices must be cleaned again.
 2. Do not re-use silicon stops.
 3. Discard devices that show any defect.
-

Packaging

- Place the devices in a kit, support or container to avoid any contact between instruments and pack the devices in sterilization pouches.
1. Avoid any contact between instruments during sterilization.
 2. For sharp devices that are not contained within a box, silicon tubes should be placed around the devices to prevent puncturing the packaging.
 3. Seal the pouches according to the recommendation of the pouch manufacturer. If a thermo-sealer is used, the process must be validated.
 4. Check the validity period of the pouch given by the pouch manufacturer to determine the shelf life.
-

Sterilization

- Sterilization at 134°C for 18 min is recommended for these devices, for the purpose of de-activating potential prions.
- The instruments must be sterilized according to the packaging labelling.
- Place the pouches in the steam sterilizer according to the recommendation given by the sterilizer manufacturer.
- Use only steam sterilizer that matches the requirements of EN 13060 (class B, small sterilizer), EN 285 (full size sterilizer).
- Use a validated sterilization procedure according to ISO 17665 with a minimum drying time of 20 min
- Respect the maintenance procedure of the sterilizer given by the sterilizer manufacturer.
- Control the efficiency and acceptance criteria of the sterilization procedure (packaging integrity, no humidity, no color change of packaging, positive physico-chemical indicators, conformity of actual cycle parameters, to reference cycle parameters).

Storage

- Keep devices in sterilization packaging in a clean environment, away from sources of moisture and direct sunlight. Store at ambient temperature.
1. Sterility cannot be guaranteed if packaging is open, damaged or wet.
 2. Check the packaging and the medical devices before using them (packaging integrity, no humidity and use by date).
 3. Storage Conditions:
 - Temperature: 10°C - 35°C Humidity: <70% RH
 - Humidity: <70% RH

8. Transportation

- To prevent damage to the medical device during transit, use a specific racks, trays or rigid containers may be recommended.
- If the package is damaged, the device is considered non-sterile and must be cleaned, disinfected, and sterilized before use.

9. Storages and duration of use

- Avoid storing at high temperature, humidity and direct sunlight. Keep away from liquids. Store at room temperature.
- Do not damage or puncture the packaging materials.
- The product is subject to be improved without previous notice. Enforce the first-in first-out method for stock management.
- Do not keep under germicidal lamp to avoid deterioration.

10. Disposal

- For proper disposal, always observe national laws and recommendations of the authorities.

11. Packaging

- Minimum packaging unit: 4pcs/in a box.
- Assortment: 1pcs from each size in one box.

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MedNet EC-REP GmbH
Wellkang Ltd
SHENZHEN ROGIN MEDICAL CO., LTD