


REVERSO SILVER protocol (treatment)



After opening the chamber and locating all canals, fill the access cavity with sodium hypochlorite.



Take a K10 file and insert it passively into the canal, rotating it 1/8th of a turn to the left and right. Move forward without applying pressure and stop when the file stops moving. Note this length.



Take your REVERSO SILVER SO25 instrument and prepare the canal to this length. To do this, insert the instrument in motion into the canal and move forward applying a very slight apical pressure. Remain in the canal for 3 seconds, then remove the file and clean the flute. Repeat this action until you have reached the desired length. You have now prepared the coronal 2/3 of the canal.



Irrigate abundantly.



Take again your K10 file and move forward passively like before. In many cases, the file will reach the working length passively. Measure this length using an apex locator. If at this point you cannot reach the working length, do not force it. Remove the K10 file and prepare the canal portion obtained using your REVERSO SILVER SO25, then repeat this step.



Prepare the canal to the working length with your REVERSO SILVER SO25 making the same movement as before. Bring your REVERSO SILVER SO25 file 3 times to the working length.

This completes the mechanical preparation of the tooth.

REVERSO SILVER protocol (retreatment)



After opening the chamber and locating all canals, fill the access cavity with sodium hypochlorite.




Use your REVERSO SILVER (or REVERSO BLUE) 25 instrument in contact with the gutta percha, applying a slight apical pressure. When the instrument engages in the gutta percha, move down a few mm and out again.



Rinse thoroughly.



Take your instrument again, put it in contact with the gutta percha and move forward a few mm.



Rinse and repeat to the working length.