

access

An **affordable** range of products
for **more cost-effective** dental practice

EASY-FILE FLEX

Variable taper


Alternating anti-clockwise
and clockwise movement



WITH HEAT
TREATMENT

The NiTi endodontic
sequence up to
50% cheaper

Test it



EASY-FILE FLEX

MONEY BACK SATISFACTION GUARANTEE*

STANDARD AND/OR
DIFFICULT CASES




WITH **EASY-FILE FLEX**


SIMPLICITY AND SPEED
ONLY 1 INSTRUMENT!




EASY-FILE FLEX




Flexible



Cutting power




Resistance



Sterilizable at 134°C

SIMPLIFICATION OF THE ENDODONTIC PROTOCOL

FOR TREATMENT




EASY-FILE FLEX	LENGTH		
	21 mm	25 mm	31 mm
ASSORTED KIT x4	205-515	205-516	205-517
S20 - 7% - ISO 20 x6	205-519	205-520	205-521
R25 - 7% - ISO 25 x6	205-523	205-524	205-525
M35 - 6% - ISO 35 x6	205-526	205-527	205-528
L45 - 5% - ISO 45 x6	205-529	205-530	205-531

Alternating anti-clockwise and clockwise movement



ALL OF THE PREPARATION WITH ONE FILE!

STANDARD AND/OR DIFFICULT CASES



* The "Money Back Satisfaction Guarantee" applies exclusively in accordance with the following terms: The guarantee is only valid within one hundred days from the date of invoice. ■ The warranty applies to unpacked and unused products, except for a box of 6 test files. ■ Products must be returned to the company in the packaging and original overpacks and accompanied by the reason for the exchange or refund request and a photocopy of a quotation or purchase invoice after having obtained an agreement to return the files to the address of the supplier. ■ In all other cases, the guarantee cannot be applied. ■ EASY-FILE FLEX: Class IIa medical device. CE 0197 marking. Notified body TÜV Rheinland. European representative: Wellkang Limited The Black Church, St. Mary's Place, Dublin 7, IRELAND. **Photos are not contractual.**