

# access

An **affordable** range of products  
for **more cost-effective** dental practice

## GO-TAPER


Variable taper

Continuous rotation



The NiTi endodontic  
sequence up to  
**50% cheaper**


**Test it**



**GO-TAPER**

**MONEY BACK SATISFACTION GUARANTEE\***

### STANDARD CASES



WITH **GO-TAPER**

**FIND  
THE ORIGINAL  
SENSATIONS**




\* The "Money Back Satisfaction Guarantee" applies exclusively in accordance with the following terms: The guarantee is only valid within one hundred days from the date of invoice. ■ The warranty applies to unpacked and unused products, except for a box of 6 test files. ■ Products must be returned to the company in the packaging and original overpacks and accompanied by the reason for the exchange or refund request and a photocopy of a quotation or purchase invoice after having obtained an agreement to return the files to the address of the supplier. ■ In all other cases, the guarantee cannot be applied.


# GO-TAPER UNIVERSAL

FIND THE ORIGINAL SENSATIONS WITH GO-TAPER


FOR TREATMENT




Flexible



Cutting power



Resistance



Sterilizable at 134°C

x6




GO-TAPER UNIVERSAL	LENGTH			
	19 mm	21 mm	25 mm	31 mm
ASSORTED KIT	-	202-268	202-276	202-284
A0 - 4% - ISO 19	202-260	-	-	-
A1 - 2% - ISO 18	-	202-261	202-269	202-277
A2 - 4% - ISO 20	-	202-262	202-270	202-278
B1 - 7% - ISO 20	-	202-263	202-271	202-279
B2 - 8% - ISO 25	-	202-264	202-272	202-280
B3 - 9% - ISO 30	-	202-265	202-273	202-281
B4 - 6% - ISO 40	-	202-266	202-274	202-282
B5 - 5% - ISO 50	-	202-267	202-275	202-283

# GO-TAPER RETREATMENT

FOR RETREATMENT

- For the **REMOVAL** of endodontic filling
- **INSTRUMENTS ADAPTED** to each third of the root canal

x6



GO-TAPER RETREATMENT	LENGTH	REFERENCE
ASSORTED KIT	Assorted	202-288
C1 - 9% - ISO 30	16 mm	202-285
C2 - 8% - ISO 25	18 mm	202-286
C3 - 7% - ISO 20	22 mm	202-287