

access

An **affordable** range of products
for **more cost-effective** dental practice

REVERSO BLUE


Variable taper

Alternating anti-clockwise
and clockwise movement



The NiTi endodontic
sequence up to
50% cheaper

Test it



REVERSO BLUE

MONEY BACK SATISFACTION GUARANTEE*

**STANDARD AND/OR
DIFFICULT CASES**

With heat treatment

FLEXIBILITY ++

**CYCLIC FATIGUE
RESISTANCE ++**

WITH **REVERSO BLUE**

**SIMPLICITY AND SPEED
ONLY 1 INSTRUMENT!**





REVERSO BLUE


WITH HEAT TREATMENT

FOR TREATMENT


x 6




Flexible ++



Cutting power



Resistance ++



Sterilizable at 134°C


REVERSO BLUE	LENGTH		
	21 mm	25 mm	31 mm
ASSORTED KIT	206-662	206-663	206-664
SO25 - 8% - ISO 25	206-666	206-668	206-669
SO40 - 6% - ISO 40	206-671	206-672	206-673
SO50 - 5% - ISO 50	206-674	206-675	206-676

Alternating anti-clockwise and clockwise movement



ALL OF THE PREPARATION
WITH ONE FILE!

STANDARD AND/OR
DIFFICULT CASES



* The "Money Back Satisfaction Guarantee" applies exclusively in accordance with the following terms: The guarantee is only valid within one hundred days from the date of invoice. ■ The warranty applies to unpacked and unused products, except for a box of 6 test files. ■ Products must be returned to the company in the packaging and original overpacks and accompanied by the reason for the exchange or refund request and a photocopy of a quotation or purchase invoice after having obtained an agreement to return the files to the address of the supplier. ■ In all other cases, the guarantee cannot be applied. ■ REVERSO BLUE: Class IIa medical device. CE 0197 marking. Notified body TÜV Rheinland. European representative: Wellkang Limited The Black Church, St. Mary's Place, Dublin 7, IRELAND. **Photos are not contractual.**