

# access

An **affordable** range of products  
for **more cost-effective** dental practice

## REVERSO SILVER


Variable taper

Alternating anti-clockwise  
and clockwise movement



The NiTi endodontic  
sequence up to  
**50% cheaper**

Test it



REVERSO SILVER

**MONEY BACK SATISFACTION GUARANTEE\***

### STANDARD CASES



WITH REVERSO SILVER

SIMPLICITY AND SPEED  
ONLY 1 INSTRUMENT!




# REVERSO SILVER


**SIMPLIFICATION OF THE  
ENDODONTIC PROTOCOL**

**FOR TREATMENT  
AND RETREATMENT**


**x 6**




Flexible



Cutting  
power



Resistance



Sterilizable  
at 134°C




REVERSO SILVER	LENGTH		
	21 mm	25 mm	31 mm
ASSORTED KIT	204-650	204-655	204-659
SO25 - 8% - ISO 25	204-651	204-656	204-660
SO40 - 6% - ISO 40	204-652	204-657	204-661
SO50 - 5% - ISO 50	204-653	204-658	204-662

Alternating anti-clockwise  
and clockwise movement



**ALL OF THE PREPARATION  
WITH ONE FILE!**

**STANDARD CASES**



\* The "Money Back Satisfaction Guarantee" applies exclusively in accordance with the following terms: The guarantee is only valid within one hundred days from the date of invoice. ■ The warranty applies to unpacked and unused products, except for a box of 6 test files. ■ Products must be returned to the company in the packaging and original overpacks and accompanied by the reason for the exchange or refund request and a photocopy of a quotation or purchase invoice after having obtained an agreement to return the files to the address of the supplier. ■ In all other cases, the guarantee cannot be applied. ■ REVERSO SILVER: Class IIa medical device. CE 0197 marking. Notified body TÜV Rheinland. European representative: Wellkang Limited The Black Church, St. Mary's Place, Dublin 7, IRELAND. **Photos are not contractual.**