

access

EASY FILL

**Gutta-percha injector
designed for obturation
of the 2/3 coronal of the canal**

- Facilitates injection of gutta-percha
- 4 preset temperatures: 100°C, 150°C, 180°C and 200°C
- Cordless gun-shaped obturation system!



Supplied with :

6 EASY FILL TIPS (2 x 25G-24MM, 2 x 25 G-28MM, 2 X 24G-24MM),
1 tip bending tool and 3 thermal insulation sleeves.

EASY FILL

3153026

EASY FILL

TIPS

134°C



Sterilisable at 134°C
for 18 minutes

x6



EASY FILL TIPS

25G (0,45MM) 24MM	3153027
25G (0,45MM) 28MM	3153028
24G (0,55MM) 24MM	3153029
23G (0,65MM) 24MM	3153030
23G (0,65MM) 28MM	3153031

access

EASY PACK

Hot condenser for root canal filling

- Used to soften gutta-percha then compact it
- 5 preset temperatures: 90°C, 150°C, 180°C, 200°C and 230°C
- Cordless small obturation system !



Supplied with:

4 assorted EASY PACK TIPS.

EASY PACK

3153032

EASY PACK

TIPS

134°C



Sterilisable at 134°C
for 18 minutes

x 4



EASY PACK TIPS

04/40	3153033
06/55	3153034
08/55	3153035
10/55	3153036