

access

STERILE R

REVERSO SILVER

Alternate movement

Clockwise and counter-clockwise




250-350 rpm


NiTi alloy

for controlled memory to prevent ledge formation and perforations

Try **REVERSO SILVER**



Simplify YOUR **ENDODONTIC PROTOCOL !**




Variable taper

STANDARD CASES



PRE - OP



Disobturation



Master cone



POST - OP




Clinical case by

Dr. Hugo Setbon

Endodontics specialist


REVERSO SILVER



Flexibility



Cutting power




Resistance



Single use

ALL THE PREPARATION WITH 1 SINGLE INSTRUMENT !

FOR TREATMENT AND RETREATMENT



x6

REVERSO SILVER	LENGTH		
	21 mm	25 mm	31 mm
ASSORTED KIT	3166602	3166606	3166610
SO25 - 8% - ISO 25	3166599	3166603	3166607
SO40 - 6% - ISO 40	3166600	3166604	3166608
SO50 - 5% - ISO 50	3166601	3166605	3166609

access

A COMPLETE RANGE OF STERILE FILES

	STANDARD CANALS STRAIGHT/WIDE	DIFFICULT CANALS NARROW/CURVED
Pre-enlargement		
	CAT-FILE	✓
	GO-GLIDER	✓
	EASY-FILE <i>Glider flexibility +++</i>	✓
Treatment		
	GO-TAPER FLEXI <i>flexibility +++</i>	✓
	GO-TAPER BLUE <i>flexibility +++</i>	✓
	GO-TAPER X	✓
	GO-TWO	✓
	GO-TAPER UNIVERSAL	✓
	REVERSO BLUE <i>flexibility +++</i>	✓
	EASY-FILE FLEX <i>flexibility +++</i>	✓
Retreatment		
	GO-TAPER RETREATMENT	✓
	REVERSO SILVER	✓