

access


STERILE R

PREFILE


FILES FOR CATHETERISATION



Try **PREFILE**



Offer yourself **MECHANISED
CATHETERISATION !**



STANDARD AND/OR
DIFFICULT CASES



- **FLEXIBILITY**
- **CYCLIC FATIGUE
RESISTANCE**
- **FOAM GUIDE TIPS**

PROTOCOL with **PREFILE**




+ preparation



PREFILE


CATHETERISATION



Flexibility



Cutting power




Resistance



Single use

ENSURE THE PERMEABILITY OF THE ROOT CANAL BEFORE ROTARY INSTRUMENTATION !



x6

PREFILE	LENGTH		
	21 mm	25 mm	31 mm
ASSORTED KIT	3148396	3148400	3148404
2% – ISO 013 STE	3148393	3148397	3148401
2% – ISO 016 STE	3148394	3148398	3148402
2% – ISO 019 STE	3148395	3148399	3148403

access

A COMPLETE RANGE OF STERILE FILES

	STANDARD CANALS STRAIGHT/WIDE	DIFFICULT CANALS NARROW/CURVED
Pre-enlargement		
	PREFILE	✓
	GO-GLIDER	✓
	EASY-FILE FLEX <i>Glider flexibility +++</i>	✓
Treatment		
	GO-TAPER FLEX <i>flexibility +++</i>	✓
	GO-TAPER BLUE <i>flexibility +++</i>	✓
	GO-TAPER X	✓
	GO-TWO	✓
	GO-TAPER UNIVERSAL	✓
	REVERSO BLUE <i>flexibility +++</i>	✓
	EASY-FILE FLEX <i>flexibility +++</i>	✓
Retreatment		
	GO-TAPER RETREATMENT	✓
	REVERSO SILVER	✓