

# access

STERILE R


## EASY-FILE FLEX *Glider*

Pre-enlarger




NiTi alloy and heat treatment

Try EASY-FILE FLEX Glider



Secure YOUR ENDODONTIC PROTOCOL !



Variable taper

STANDARD CASES



FLEXIBILITY



ENDURANCE



FOAM GUIDE TIPS

POST - OP


Clinical case by  
**Dr. Hugo Setbon**  
Endodontics specialist

with **EASY-FILE FLEX Glider**  
then EASY-FILE FLEX  
for the preparation

PRE - OP




# EASY-FILE FLEX Glider



Flexibility



Cutting power




Resistance



Single use

TO ANTICIPATE A CONTROLLED ROOT CANAL PREPARATION !

PRE-ENLARGER



x6

EASY-FILE FLEX Glider	LENGTH		
	21 mm	25 mm	31 mm
2% – ISO 015	3148405	3148406	3148407

access

A COMPLETE RANGE OF STERILE FILES

		STANDARD CANALS STRAIGHT/WIDE	DIFFICULT CANALS NARROW/CURVED
<b>Pre-enlargement</b>			
	<b>PREFILE</b>		✓
	<b>GO-GLIDER</b>		✓
	<b>EASY-FILE FLEX Glider</b> flexibility +++		✓
<b>Treatment</b>			
	<b>GO-TAPER FLEX</b> flexibility +++	✓	✓
	<b>GO-TAPER BLUE</b> flexibility +++	✓	✓
	<b>GO-TAPER X</b>	✓	✓
	<b>GO-TWO</b>	✓	✓
	<b>GO-TAPER UNIVERSAL</b>	✓	✓
	<b>REVERSO BLUE</b> flexibility +++	✓	✓
	<b>EASY-FILE FLEX</b> flexibility +++	✓	✓
<b>Retreatment</b>			
	<b>GO-TAPER RETREATMENT</b>	✓	
	<b>REVERSO SILVER</b>	✓	