

access

STERILE R

EASY-FILE FLEX

Variable taper

STANDARD AND/OR
DIFFICULT CASES

Alternate movement

Clockwise and counter-clockwise



250-350 rpm

NiTi alloy and heat treatment

for controlled memory to prevent
ledge formation and perforations

Try EASY-FILE FLEX




Simplify YOUR ENDODONTIC PROTOCOL
WITH A **SINGLE INSTRUMENT** !



PRE - OP



Master cone




POST - OP



Clinical case by
Dr. Hugo Setbon
Endodontics specialist


EASY-FILE FLEX



Flexibility



Cutting power




Resistance



Single use

ALL THE PREPARATION WITH 1 SINGLE INSTRUMENT !

FOR TREATMENT



EASY-FILE FLEX	LENGTH		
	21 mm	25 mm	31 mm
ASSORTED KIT x4	3166615	3166620	3166625
S20 - 7% - ISO 20 x6	3166611	3166616	3166621
R25 - 7% - ISO 25 x6	3166612	3166617	3166622
M35 - 6% - ISO 35 x6	3166613	3166618	3166623
L45 - 5% - ISO 45 x6	3166614	3166619	3166624

access

A COMPLETE RANGE OF STERILE FILES

	STANDARD CANALS STRAIGHT/WIDE	DIFFICULT CANALS NARROW/CURVED
Pre-enlargement		
	CAT-FILE	✓
	GO-GLIDER	✓
	EASY-FILE Glider flexibility +++	✓
Treatment		
	GO-TAPER FLEXI flexibility +++	✓
	GO-TAPER BLUE flexibility +++	✓
	GO-TAPER X	✓
	GO-TWO	✓
	GO-TAPER UNIVERSAL	✓
	REVERSO BLUE flexibility +++	✓
	EASY-FILE FLEX flexibility +++	✓
Retreatment		
	GO-TAPER RETREATMENT	✓
	REVERSO SILVER	✓