

# access

STERILE R

## GO-GLIDER

Pre-enlarger




NiTi alloy

for optimal fracture resistance

Try GO-GLIDER



Regain FLEXIBILITY AND ENDURANCE !



Variable taper

STANDARD CASES



- FLEXIBILITY
- ENDURANCE
- FOAM GUIDE TIPS


### PROTOCOL with GO-GLIDER



+ preparation




# GO-GLIDER



Flexibility



Cutting power




Resistance



Single use

TO ANTICIPATE A CONTROLLED ROOT CANAL PREPARATION !

PRE-ENLARGER



x6

GO-GLIDER	LENGTH		
	21 mm	25 mm	31 mm
2% - ISO 016	3166663	3166664	3166665

## access

A COMPLETE RANGE OF STERILE FILES

	STANDARD CANALS STRAIGHT/WIDE	DIFFICULT CANALS NARROW/CURVED
<b>Pre-enlargement</b>		
	<b>CAT-FILE</b>	✓
	<b>GO-GLIDER</b>	✓
	<b>EASY-FILE</b> <i>Glider flexibility +++</i>	✓
<b>Treatment</b>		
	<b>GO-TAPER FLEXI</b> <i>flexibility +++</i>	✓
	<b>GO-TAPER BLUE</b> <i>flexibility +++</i>	✓
	<b>GO-TAPER X</b>	✓
	<b>GO-TWO</b>	✓
	<b>GO-TAPER UNIVERSAL</b>	✓
	<b>REVERSO BLUE</b> <i>flexibility +++</i>	✓
	<b>EASY-FILE FLEX</b> <i>flexibility +++</i>	✓
<b>Retreatment</b>		
	<b>GO-TAPER RETREATMENT</b>	✓
	<b>REVERSO SILVER</b>	✓