

# access

STERILE R

## GO-TAPER UNIVERSAL

Continuous rotation



250-300 rpm

NiTi alloy

for optimal fracture resistance

Try GO-TAPER




Optimize DEBRIS REMOVAL !



Variable taper

STANDARD CASES



PRE - OP




POST - OP



Clinical case by  
**Dr. Hugo Setbon**  
Endodontics specialist


# GO-TAPER UNIVERSAL



Flexibility



Cutting power



Resistance




Single use

**SAME SEQUENCE FOR ALL STANDARD CANALS !**

FOR TREATMENT

x6



GO-TAPER UNIVERSAL	LENGTH			
	19 mm	21 mm	25 mm	31 mm
ASSORTED KIT	-	3166578	3166586	3166594
A0 - 4% - ISO 19	3166570	-	-	-
A1 - 2% - ISO 18	-	3166571	3166579	3166587
A2 - 4% - ISO 20	-	3166572	3166580	3166588
B1 - 7% - ISO 20	-	3166573	3166581	3166589
B2 - 8% - ISO 25	-	3166574	3166582	3166590
B3 - 9% - ISO 30	-	3166575	3166583	3166591
B4 - 6% - ISO 40	-	3166576	3166584	3166592
B5 - 5% - ISO 50	-	3166577	3166585	3166593

## access

A COMPLETE RANGE OF STERILE FILES

	STANDARD CANALS STRAIGHT/WIDE	DIFFICULT CANALS NARROW/CURVED
<b>Pre-enlargement</b>		
	<b>CAT-FILE</b>	✓
	<b>GO-GLIDER</b>	✓
	<b>EASY-FILE</b> <i>Glider flexibility +++</i>	✓
<b>Treatment</b>		
	<b>GO-TAPER FLEXI</b> <i>flexibility +++</i>	✓
	<b>GO-TAPER BLUE</b> <i>flexibility +++</i>	✓
	<b>GO-TAPER X</b>	✓
	<b>GO-TWO</b>	✓
	<b>GO-TAPER UNIVERSAL</b>	✓
	<b>REVERSO BLUE</b> <i>flexibility +++</i>	✓
	<b>EASY-FILE FLEX</b> <i>flexibility +++</i>	✓
<b>Retreatment</b>		
	<b>GO-TAPER RETREATMENT</b>	✓
	<b>REVERSO SILVER</b>	✓