

access

STERILE R

GO-TAPER BLUE

Continuous rotation




250 - 350 rpm


NiTi alloy and heat treatment

for controlled memory to prevent
ledge formation and perforations

Try GO-TAPER BLUE



Regain FLEXIBILITY AND ENDURANCE !



Variable taper

STANDARD AND/OR
DIFFICULT CASES



PRE - OP




POST - OP



Clinical case by
Dr. Hugo Setbon
Endodontics specialist


GO-TAPER BLUE



Flexibility



Cutting power



Resistance




Single use

NEW THERMAL TREATMENT "BLUE TREATMENT"

x6

FOR TREATMENT



GO-TAPER BLUE	LENGTH			
	19 mm	21 mm	25 mm	31 mm
ASSORTED KIT A0 à B3	–	3166634	3166642	3166650
A0 – 4% – ISO 019	3166626	–	–	–
A1 – 2% – ISO 018	–	3166627	3166635	3166643
A2 – 4% – ISO 020	–	3166628	3166636	3166644
B1 – 7% – ISO 020	–	3166629	3166637	3166645
B2 – 8% – ISO 025	–	3166630	3166638	3166646
B3 – 9% – ISO 030	–	3166631	3166639	3166647
B4 – 6% – ISO 040	–	3166632	3166640	3166648
B5 – 5% – ISO 050	–	3166633	3166641	3166649

access

A COMPLETE RANGE OF STERILE FILES

	STANDARD CANALS STRAIGHT/WIDE	DIFFICULT CANALS NARROW/CURVED
Pre-enlargement		
	CAT-FILE	✓
	GO-GLIDER	✓
	EASY-FILE <i>Glider flexibility +++</i>	✓
Treatment		
	GO-TAPER FLEXI <i>flexibility +++</i>	✓
	GO-TAPER BLUE <i>flexibility +++</i>	✓
	GO-TAPER X	✓
	GO-TWO	✓
	GO-TAPER UNIVERSAL	✓
	REVERSO BLUE <i>flexibility +++</i>	✓
	EASY-FILE FLEX <i>flexibility +++</i>	✓
Retreatment		
	GO-TAPER RETREATMENT	✓
	REVERSO SILVER	✓