

access

STERILE R


GO-TAPER X

Continuous rotation 
250-300 rpm

NiTi alloy

for controlled memory to prevent
ledge formation and perforations

Try GO-TAPER X



Optimize DEBRIS REMOVAL !



Variable taper

STANDARD CASES

CLINICAL CASE 1

PRE - OP



POST - OP



CLINICAL CASE 2

PRE - OP




POST - OP



Clinical case by
Dr. Hugo Setbon
Endodontics specialist


GO-TAPER X



Flexibility



Cutting power



Resistance




Single use

FELXIBILITY AND ENDURANCE !

Off-centre rectangular cross-section design

Asymmetrical rotary movement

FOR TREATMENT



x6

GO-TAPER X	LENGTH		
	21 mm	25 mm	31 mm
ASSORTED KIT	3166711	3166717	3166723
4% – ISO 017	3166706	3166712	3166718
6% – ISO 017	3166707	3166713	3166719
7% – ISO 030	3166708	3166714	3166720
6% – ISO 040	3166709	3166715	3166721
6% – ISO 050	3166710	3166716	3166722

access

A COMPLETE RANGE OF STERILE FILES

	STANDARD CANALS STRAIGHT/WIDE	DIFFICULT CANALS NARROW/CURVED
Pre-enlargement		
	CAT-FILE	✓
	GO-GLIDER	✓
	EASY-FILE <i>Glider flexibility +++</i>	✓
Treatment		
	GO-TAPER FLEXI <i>flexibility +++</i>	✓
	GO-TAPER BLUE <i>flexibility +++</i>	✓
	GO-TAPER X	✓
	GO-TWO	✓
	GO-TAPER UNIVERSAL	✓
	REVERSO BLUE <i>flexibility +++</i>	✓
	EASY-FILE FLEX <i>flexibility +++</i>	✓
Retreatment		
	GO-TAPER RETREATMENT	✓
	REVERSO SILVER	✓